

# A Breathable Barrier Face Mask by Hudson Overall Company, Inc. in Partnership with GD Apparel Made with Burlington® Maxima® Fabric



**Our masks use fabrics that protect.**

Burlington® Performance Fabrics is committed to engineering and producing state-of-the-art fabrics that combine the highest degree of protection, performance, comfort, and durability.

** Hudson Overall Company, Inc.  GD Apparel  Burlington®**

Hudson Overall Company, Inc. is Dedicated to Building #GoodsWithPurpose  
using Quality Materials & Together with GD Apparel we are Practicing Manufacturing Excellence.

Maxima® fabrics should only be used for minimum and low performance category face masks. These fabrics should not be used in the production of N95 masks. When used over top of an N95 masks you may be able to increase the life of the N95 masks and reuse it per CDC guidelines. Neither Burlington, Hudson Overall Co., or GD Apparel makes representations or warranties about the level of protection from viral or bacterial pathogens. These products are intended for personal care and industrial use only, and the aforementioned parties make no claim that their use can reduce or prevent the transmission of disease or illness.

## Breathable Barrier Fabric We Use

- Maxima® AT – style 4041 – Medical fabric constructed with non-linting polyester fibers and fluid repellent finish made from air textured yarns for natural hand and better comfort.
- Maxima® AT ESD – style 4045 – Anti-static version of Maxima® AT, made from air textured yarns with a carbon stripe and fluid repellent finish for a natural hand and great comfort.

### Fabric: Liquid Repellent Finish Breathable Barrier Fabric



**Trim:** 3/16" 3-cord latex-free elastic    **Thread:** Tex-40

- Single-layer design permits breathability and air filtration during respiration utilizing high-tech barrier fabric from Burlington Fabrics, and therefore does not require a liner for improvement of breathability. You breathe through this mask, not around it like other product offerings available in the market.
- The MAXIMA® AT barrier finish treatment lifespan lasts up to 75 sterile laundry cycles (~160°F). Product lifespan can be elongated with correct hand-washing maintenance.
  - To machine wash, simply place through one hot wash cycle, and hang dry or pat dry with clean cloth or towel.
  - To sink wash, use hot water, and apply anti-bacterial soap on both sides, covering liberally. Scrub well, then rinse.
  - Pat dry with paper towel or cloth or use wall unit to dry. Note: fabric is liquid intolerant, and will dry quickly, therefore elastic is the main focus of drying process.
- Three sizes available\*: S/M (red stitching) L/XL (blue stitching), 2XL/3XL (green stitching)  
\*NOTE: For use over N95 respirator, we recommend sizing up one full size.

Burlington® and Maxima® are registered trademarks of Burlington Industries LLC, a subsidiary of Elevate Textiles, Inc.

Maxima® fabrics should only be used for minimum and low performance category face masks. These fabrics should not be used in the production of N95 masks. When used over top of an N95 masks you may be able to increase the life of the N95 masks and reuse it per CDC guidelines. Neither Burlington, Hudson Overall Co., or GD Apparel makes representations or warranties about the level of protection from viral or bacterial pathogens. These products are intended for personal care and industrial use only, and the aforementioned parties make no claim that their use can reduce or prevent the transmission of disease or illness.